

AutoSave Off | HLA XLS in NMDP Code + Ambiguous Combos + Common Combo | Search

File Home Insert Draw Page Layout Formulas Data Review View Help | Share | Comments

SampleID	A1	A2	B1	B2	C1	C2	DRB11	DRB12	DRB31	DRB32	DRB41	DRB42	DRB51	DRB52	DQB11	DQB12	DQA11	DQA12	DPB11	DPB12	DPA11	DPA12	HLA-A	HLA-B	HLA-C	HLA-DRB1
Sample ID 2	'0000	'0000	'0000	'0000	'0000	'0000	'0000	'0000							'0000	'0000			'0000	'0000			A*0000/A*0000	B*0000/B*0000	C*0000/C*0000	DRB1*0000/DRB1*0000

Report columns will be continued for each sample as shown below:

SampleID	HLA-DRB3	HLA-DRB4	HLA-DRB5	HLA-DQB1	HLA-DQA1	HLA-DPB1	HLA-DPA1	HLA-A	HLA-B	HLA-C	HLA-DRB1	HLA-DRB3	HLA-DRB4	HLA-DRB5	HLA-DQB1	HLA-DQA1	HLA-DPB1	HLA-DPA1	Comments
Sample ID 2				DQB1*0000/DQB1*0000			DPB1*0000/DPB1*0000												

ARS - Antigen Recognition Site
 Reference - Human Immunology 68, 392-417 (2007) @American Society for Histocompatibility and Immunogenetics, 2007

Note:
 Alleles bearing suffix 'G' in Class I locus has identical sequences in exon 2 and exon 3 antigen recognition sites.
 Alleles bearing suffix 'G' in Class II locus has identical sequences in exon 2 antigen recognition sites.
 Please find below the list of alleles included for each allele suffix with 'G'.
 Allele Database Version used in the report is 3.34.0(Oct 2018)
 Sequencing Technology : Illumina Sequencing
 NMDP code is reported if G code is not available after analyzing non-ARS exons

ActualTyping | Appendix

AutoSave Off | HLA XLS in NMDP Code + Ambiguous Combos + Common Combo | Search

File Home Insert Draw Page Layout Formulas Data Review View Help | Share | Comments

Sample ID	Locus	NMDP code	G code	Segment Sequenced	Allele Database Version	Included Alleles

ActualTyping | Appendix